

Sievers* M500 Total Organic Carbon Analyzer

Accuracy. Efficiency. Integrity.



Advancing the Sievers Legacy

The Sievers M500 Online Total Organic Carbon (TOC) Analyzer revolutionizes online detection of organics by bringing new performance, design, and data management features to the industry-leading Sievers TOC platform.

Advancing the Sievers Legacy

The Sievers reagentless Membrane Conductometric TOC Detection Method delivers unmatched reliability, robustness, and accuracy for pharmaceutical Purified Water (PW) and Water for Injection (WFI) applications.

Sievers products have been leaders in TOC analytical instrumentation since the introduction of membrane conductometric technology to the market more than 25 years ago.

Superior Performance and Features

The Sievers M500 is our third generation online TOC analyzer designed for accuracy, efficiency, and integrity.

Building upon the proven results of the 500 RL, the M500 improves our industry-leading performance and adds a suite of cutting-edge features:

- Super iOS 4-port sampler automates grab samples to increase efficiency.
- Reducing analysis time by 50% facilitates utilization of real-time data, early detection, and process control.
- 10-inch touchscreen enables faster and easier setup and operation.
- Standardized and customized protocols improve productivity.
- Enhanced data features improve data integrity.

50%

Reduction in
Analysis Time

Designed for: Data Integrity

Sievers M500: Designed for Today's Data

The Sievers M500 leverages comprehensive industry-leading TOC data management tools to ensure data is secure.

The Sievers M500 supports compliance to 21 CFR Part 11 and adherence to US FDA and other Pharmacopoeia Data Integrity guidelines using new digital features:

Data Transfer

- Remote access
- Ethernet and WiFi
- Advanced communications using 4-20 mA, Modbus, Profinet, and binary

Data Security


- Password protection
- Data encryption
- Customizable access, roles, and permissions

Web-Based Data Management

- Closed system architecture
- Customizable data transfer and export

Sievers Membrane Conductometric TOC Detection Technology

Proprietary gas-permeable membrane selectively passes only the CO_2 produced from oxidized organics



Unmatched reliability, robustness, and accuracy


Higher selectivity, sensitivity, stability, and precision

Discrimination between **organic** and **inorganic** carbon provides critical information for process control

Legend


- Organic Carbon Molecule
- SO_4^{2-}
- Cl^-
- CO_2
- HCO_3^-
- H^+

Sievers M500



Accuracy

- Range: 0.03 ppb – 2.5 ppm
- Limit of Detection: 0.03 ppb
- Limit of Quantification: 0.1 ppb
- Higher selectivity prevents interference from other compounds
- Measurements compliant to USP, EP, JP, IP, ChP, KP and all other harmonized compendia



Efficiency

- 50% reduction in analysis time
- Simultaneous TOC & conductivity measurement
- Super iOS 4-port sampler
- Automated system suitability, calibration, & validation
- Adaptive auto-zero



Integrity

- Designed for 21 CFR Part 11
- Remote access
- Advanced communications with Ethernet & WiFi
- Improved data transfer, data security, and data management



Key Features

Compliant

- US Pharmacopeia (USP)
- European Pharmacopoeia (EP)
- Japanese Pharmacopoeia (JP)
- Indian Pharmacopoeia (IP)
- Chinese Pharmacopoeia (ChP)
- Korean Pharmacopoeia (KP)
- All other harmonized compendia for TOC and conductivity analysis

Super iOS

- 4-port sampler
- Automated calibration and protocols
- 2x faster results
- 50% faster validation protocols
- Reduced human error and labor costs

10-inch touchscreen

- Intuitive interface
- At-a-glance viewing
- Faster and easier setup & operation
- Chemical resistant



Data Integrity

- Supports 21 CFR Part 11 & Data Integrity regulations and guidance
- Lean compendia water analysis
- Simultaneous TOC & conductivity measurement

Advanced communications

- Ethernet, WiFi, 4-20 mA, Modbus, Profinet, Profibus, and Serial connectivity

For more information, visit watertechnologies.com/sievers



Full Lifecycle Support and Services

Veolia provides a complete lifecycle solution for Sievers products. We offer a suite of tools and services to optimize your instrument from installation and validation, to operation and maintenance.

Trust Sievers validation services and Preventive Maintenance Agreements to provide peace of mind by ensuring optimal performance and providing the documentation required to support compliance, quality, and regulatory needs.

Our support team provides remote and onsite support for Sievers instruments. We are uniquely positioned to offer Failure Analysis Reports (FARs) to help close Out-of-Specification investigations faster and more efficiently.



Installation

- Field Service Engineers (FSE)
- Onsite installation
- Training
- Documentation
- Software

Validation

- Field Service Engineers (FSE)
- Validation Packages
- Onsite validation services

Operation

- Certified reference materials
- Vials and consumables
- Customized reference standards
- Preventive Maintenance Agreements (PMA)

Support

- Technical support
- Diagnostics
- Onsite repairs
- Failure Analysis Reports (FAR)
- Warranty
- Upgrades



Ask about Sievers Certified Plus full coverage service and repair

Standards

Accredited Certified Reference Materials
Standard Protocols
Low TOC Vials



Customized Reference Standards
Customized Protocols
Dual Use Conductivity & TOC (DUCT) Vials

Customized Solutions



Certified Reference Materials and Customized Solutions

In addition to providing a wide range of accredited certified reference materials, Veolia provides custom reference standards for your unique applications.

The Sievers M500 includes an expanded portfolio of instrument protocols as well as options to customize protocols for your application.

Specifications

| TOTAL ORGANIC CARBON | |
|-------------------------------------|---|
| Linear Range | 0.03 to 2,500 ppb as TOC |
| Accuracy | ± 5% of measurement; ±0.1 ppb |
| Precision | ± 1% of measurement; ± 0.03 ppb |
| Analysis Modes | Online, Online Averaged, Online Timed, Grab |
| Analysis Time ^{1,2} | 3 minutes for continuous online measurements |
| Ozone Compatibility | 50 ppb O ₃ continuous; 200 ppb O ₃ for 2 hours daily |
| Sample Flow Rate (nominal) | 0.25 mL/min |
| External Flow Rate | Minimum 50 mL/min |
| Sample Temperature | 1-95°C (34-203°F) (withstands short-term steam exposure) |
| Sample Pressure | Up to 100 psig |
| Interferences | Insensitive to organic heteroatoms |
| Calibration Stability | Typically stable for 12 months |
| Display Readout | 3 significant digits |
| CONDUCTIVITY | |
| Range | 0.01 to 800 µS/cm |
| Accuracy | 0.005 µS/cm or 1%, whichever is larger |
| Precision | ≤ 1.0 % RSD |
| INSTRUMENT | |
| Power Requirements | 100–240 VAC, 70 W, 50/60 Hz |
| Fuses | No user-replaceable fuses |
| Normal Operating Environment | Intended for indoor use only |
| Ambient Temperature | 5-40°C (41-104°F) |
| Maximum Relative Humidity | Up to 95%, noncondensing |
| Maximum Altitude | 3,000 m (9,843 ft) |
| Inputs | Two isolated binary inputs |
| Outputs | Serial (RS-232), one USB, three 4-20 mA, four alarms, one Ethernet |
| Installation/Overvoltage Category | II (protects against transients present in Category II power) |
| Safety Certifications | CE, ETL listed. Conforms to UL Std. 61010-1. Certified to CSA 22.2 No. 61010-1. |
| Pollution Degree | 2 (normally only non-conductive pollution) |
| Display | Backlit 10.1 in, 1280 x 800, touchscreen display |
| Size | H: 43.4 cm (17.1 in) W: 55.9 cm (22.0 in) D: 28.7 cm (11.3 in) |
| Weight | 16.3 kg (36 lb) |
| IP Rating | IP 55 |
| Optional Wi-Fi ³ | 802.11ac/a/b/g/n Dual band 2.4/5 GHz |
| Industrial Communications Protocols | Modbus TCP/IP & Profinet (Included), Profibus (Optional) |

¹Time to first measurement is 10 minutes.

² Base model analysis time is six minutes for continuous online measurement.

³ Not available in all countries

This information herein may be subject to change without notice and is provided for general guidance only. The dimensions and performance of systems, products and services may vary. Pictures are for example purposes and not to scale. All legal obligations are exclusively as set out in contractual documents. Nothing contained herein constitutes a representation, warranty or undertaking.

6060 Spine Road
 Boulder, CO 80301-3687 USA
 T +1 800 255 6964
 T +1 303 444 2009
 F +1 303 527 1797



Find a contact near you by visiting watertechnologies.com/sievers and clicking on “Contact Us”.

www.watertechnologies.com/sievers

